

Medical Supplies

Monitor

Monitor's
- Accessories and Consumables -

Monitor Disposable Electrodes

Vitrode *The result of over 30 years experience in disposable electrode manufacturing*

Disposable electrodes, Vitrode L

G207 (L-150X)



- 35 mm diameter for general adult use
- For general ECG and 12-lead ECG monitoring
- No gel remains when removed from skin
- Low impedance for clear waveforms
- Gentle on the skin
- Radiolucent

Disposable electrodes, Vitrode F

G210D (F-150M)




G210C (F-150S)



- Smaller size for neonate to child
- Entire attachment area has conductive and adhesive gel to reduce skin irritation
- Contain a bitter agent to help prevent accidental injection by children

Disposable electrodes, Vitrode V

G272C (V-060M6)




- Smaller size for neonate to child
- Lead wires are attached together to prevent tangling in incubators

Code (Model)	Qty	Usage	Size	Cord
G207 (L-150X)	30 pcs/pack × 5	Adult, radiolucent	35 mm dia	–
G210D (F-150M)	3 pcs/pack × 50	Adult, whole surface adhesion	25 × 45 mm	–
G210C (F-150S)	3 pcs/pack × 50	Neonate to child, whole surface adhesion	18 × 36 mm	–
G272C (V-060M6)	6 colors/pack × 10	Adult/child, prewired	25 × 45 mm	1.0 m


Disposable electrodes, Vitrode N

G300D (N-01IS3)



Designed for High-Humidity Incubators

- Durable and adhesive gel hardly swells, preventing electrodes from peeling even under humid environment
- Smallest electrodes among Vitrode series
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby



Adhesion		Skin Friendliness	Adhesive Strength to Body Movement
Humidified Environment*	Non-Humidified Environment		
★★★★★	★★★	★★★★★	★★

G300A (N-03IS3)



Skin-Friendly with Gentle Adhesive

- Skin-friendly adhesive gel; gentle for less skin irritation
- Electrodes for the baby whose dermis is immature
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby

Adhesion		Skin Friendliness	Adhesive Strength to Body Movement
Humidified Environment*	Non-Humidified Environment		
★	★★	★★★★★	★

Disposable electrodes, Vitrode V

G271A (V-120S3)



Available from Neonate to Infant

- Bitter chemical compound is used on the electrodes to prevent accidental swallowing
- Electrodes are specifically designed for the more active baby
- Soft radiolucent lead wires are bundled to prevent the leads being tangled around the baby

Adhesion		Skin Friendliness	Adhesive Strength to Body Movement
Humidified Environment*	Non-Humidified Environment		
★★	★★★★★	★★	★★★★★

*Humidified Environment means inside incubator; indicates durability of electrodes in incubator.

• Star indicates the strength in each category. The more colored stars (★) the better the electrodes address the category.

Code (Model)	Qty	Usage	Size	Cord
G271A (V-120S3)	3 colors/pack × 40	Infant to neonate, prewired	20 × 20 mm	0.6 m
G300A (N-03IS3)	3 colors/pack × 10	NICU, prewired	14 × 25 mm	0.6 m
G300D (N-01IS3)	3 colors/pack × 10	Very low birth weight (less than 1500 g) Extremely low birth weight (less than 1,000 g)	15 × 15 mm	0.6 m

Monitor ECG Cables, Leads

Electrode leads



for 3 electrodes

Code (Model)	Type	Lead
K911 (BR-903P)	IEC	Clip
K911A (BR-903PA)	AHA	Clip
K910A (BR-913P)	IEC	Snap
K910B (BR-913PA)	AHA	Snap



for 6 electrodes

Code (Model)	Type	Lead
K912 (BR-906P)	IEC	Clip
K912A (BR-906PA)	AHA	Clip
K915 (BR-916P)	IEC	Snap
K915A (BR-916PA)	AHA	Snap



for 2 electrodes

Code (Model)	Lead
K907A (BR-902P)	Clip
K908A (BR-912P)	Snap



for 3 electrodes

Code (Model)	Type	Lead
K916 (BR-963P)	IEC	Clip

ECG connection cords



for BSM/CSM/PVM, 3/6 electrodes

Code (Model)	Type	Cord
K922 (JC-906P)	IEC	3 m
K925 (JC-916P)	IEC	1.5 m
K922A (JC-906PA)	AHA	3 m

ECG patient cables



for BSM/CSM, 10 electrodes

Code (Model)	Type	Cord
K901 (BJ-900P)	IEC	3.8 m
K901A (BJ-900PA)	AHA	3.8 m

ECG connection cords



for BSM/CSM, 10 electrodes

Code (Model)	Type	Cord
K921 (JC-900P)	IEC	3 m
K921A (JC-900PA)	AHA	3 m

*Use with Vitrode V











NEONATAL adapter



Code (Model)
L021 (JC-990P)

BLUPRO *Reliable SpO₂ probes from the company that invented SpO₂*

Usage

		Low birth weight neonate	Neonate	Infant	Child			Adult		
		[Weight]	1000g	3kg	10kg	20kg	30kg	40kg	50kg	
One-patient use	TL-271T 							[Orange bar from 40kg to 50kg]		
	TL-272T 						[Orange bar from 30kg to 50kg]			
	TL-273T 			[Orange bar from 3kg to 10kg]					[Orange bar from 40kg to 50kg]	
	TL-274T 				[Orange bar from 10kg to 20kg]					
	TL-260T 	Tape S instep	Tape L instep		Tape S finger/toe			[Orange bar from 40kg to 50kg with P256 ear clip]		
	TL-535U 	Tape S	Tape L							
Reusable	TL-201T 						[Orange bar from 30kg to 50kg]			
	TL-631T3 						[Orange bar from 30kg to 50kg]			
	TL-630T3 								[Orange bar from 50kg to 50kg]	
	TL-220T 			[Orange bar from 3kg to 50kg]						


Nihon Kohden recommends combination use of one-patient use and reusable probes depending on patient condition.

BLUPRO SpO₂ probes

from the company that invented SpO₂

One-Patient use SpO₂ probes

Self-adhesive type




Adult finger/toe

Code (Model)	Cord	Qty	Code (Model)	Cord	Qty
P203A (TL-271T)	0.8 m	24	P204A (TL-271T)	0.8 m	5
P203E (TL-271T3)	1.6 m	24	P204E (TL-271T3)	1.6 m	5

Child finger/toe

Code (Model)	Cord	Qty	Code (Model)	Cord	Qty
P203B (TL-272T)	0.8 m	24	P204B (TL-272T)	0.8 m	5
P203F (TL-272T3)	1.6 m	24	P204F (TL-272T3)	1.6 m	5



Neonate instep

Code (Model)	Cord	Qty	Code (Model)	Cord	Qty
P203C (TL-273T)	0.8 m	24	P204C (TL-273T)	0.8 m	5
P203G (TL-273T3)	1.6 m	24	P204G (TL-273T3)	1.6 m	5

Infant finger/toe

Code (Model)	Cord	Qty	Code (Model)	Cord	Qty
P203D (TL-274T)	0.8 m	24	P204D (TL-274T)	0.8 m	5
P203H (TL-274T3)	1.6 m	24	P204H (TL-274T3)	1.6 m	5



Non-adhesive, for sensitive skin

Code (Model)	Patient	Qty
P205A (TL-260T)	Neonate Adult (with ear clip)	1 sensor + 4 sponge tape L/set x 5



Attachment tape for TL-260T

Code (Model)	Size	Qty
P260E (-)	S: 1 kg or less at instep 3 kg to 10 kg at finger/toe	24
P260F (-)	L: 1 kg or more	24

Ear clip, adult use with TL-260T


Code (Model)
P256 (-)

Non-adhesive, for sensitive skin

Code (Model)	Attachment site	Qty
P206 (TL-535U)	Neonates and preterm infant	1 sensor + 3 sponge tape S/set x 5

Attachment tape for TL-535U

Code (Model)	Size	Qty
P264A (-)	S: 1 kg or less	24
P264B (-)	L: 2.5 kg or less	24



BluPRO SpO₂ probes

from the company that invented SpO₂

Reusable SpO₂ probes



Finger-clip type, 20 kg or more

Code (Model)	Cord
P225F (TL-201T)	1.6 m
P225H (TL-201T)	0.6 m



Multi-site, from neonate to adult, finger or instep

Code (Model)
P225G (TL-220T)



Finger-tip type, 20 kg or more

Code (Model)
P311C (TL-631T3)

Finger-tip type, 50 kg or more

Code (Model)
P310C (TL-630T3)

Fixing tapes



Code (Model)		Qty
P259 (-)	COTTONY tape	20
P263 (-)	BluPRO attachment tape	90
P267 (-)	Probe fastener	30

SpO₂ adapters

SpO₂ adapters for P206 (TL-535U)



for BSM/CSM/PVM

Code (Model)
K931 (JL-900P)



for OLV

Code (Model)
L311 (JL-302T)



with BSM/CSM/PVM

Code (Model)
K930 (JL-930P)



with K931 (JL-900P) or L311 (JL-302T)

Code (Model)	Lead
K929A (JL-030U1)	0.3 m
K929B (JL-030U2)	1 m

Monitor CO₂


CO₂ Sensor



*Nihon Kohden's unique sensing technology
both intubated and non-intubated mainstream CO₂ measurement*



- Weighs only 4 grams
- IPX 7 protection against immersion in water
- Accurate for intubated and non-intubated patients
- Reliable even under high humidity conditions
- Durable enough for use in harsh environments
- Safe heater-less system prevents burn injuries



Code (Model)
P910A (TG-980P)



Airway Adapter for intubated patients

R805 (YG-211T)



R806 (YG-213T)

Code (Model)	Dead Space	Weight	Usage	Connection	Qty
R805 (YG-211T)	4 mL	>7 kg	Adult	-	30
R806 (YG-213T)	0.5 mL	-	Infant/ neonate	Connecting with flow sensor or respiration circuit without inner connector	30
R807 (YG-214T)	1.8 mL	2.5 to 7 kg	Infant/ neonate	Connecting flow sensor with an inner connector	30

cap-ONE mask for non-intubated patients



V938A (YG-272T)




V938C (YG-282T)



V933C (YG-232T)



V935C (YG-242T)



Code (Model)	Dead Space	Weight	Usage	Qty
V938A (YG-272T)	7 mL	30 kg or more	Adult	10
V938C (YG-282T)	10 mL	40 kg or more	Large adult	10
V933C (YG-232T)	3.5 mL	20 - 40 kg	Pediatric	10
V935C (YG-242T)	2.5 mL	7 - 20 kg	Infant	10

Ensure Quality of care under Sedation


- Accurate respiratory monitoring with optimal oxygen supply

Biteblock during Endoscopy



Code (Model)	Usage	Qty
V939A (YG-227T)	Adult	30

Alarm response time



CO₂ measurement detects respiration failure earlier than SpO₂ and ECG measurement. CO₂ monitoring leads early intervention for patient safety.

CO₂ selection guide

		P907 (TG-920P)	P903 (TG-900P)	P910A (TG-980P)	GF-210R	GF-220R
Multigas		—	—	—	✓	✓
Flow/Paw		—	—	—	—	✓
Sidestream		—	—	—	✓	✓
Mainstream	Semi-quantitative	✓	✓	—	—	—
	Quantitative	—	—	✓		
Intubated		✓	✓	✓	✓	✓
Non-intubated		✓	—	✓	—	—
Neonate		—	—	✓	—	—

✓: Available —: Not available

CO₂ mainstream, intubated and non-intubated (semi-quantitative method)

CO₂ sensor kit




Code (Model)
P907* (TG-920P)
*TG-121T + JG-920P

CO₂ adapter



Code (Model)
K984 (JG-920P)

Nasal adapters (non-intubated)



Code (Model)	Dead Space	Weight	Qty
V922 (YG-121T)	1.2 mL	>10 kg	30
V923* (YG-122T)	1.2 mL	>10 kg	30

*with O₂ supply tube adapter

Airway adapter (intubated)



Code (Model)	Dead Space	Weight	Qty
R804 (YG-111T)	4 mL	>7 kg	30

CO₂ mainstream, intubated (semi-quantitative method)

CO₂ sensor kit



Code (Model)
P903* (TG-900P)
*TG-101T + JG-900P

CO₂ adapter



Code (Model)
K981 (JG-900P)

Airway adapter




Code (Model)	Dead Space	Weight	Qty
R801 (YG-101T)	5 mL	>10 kg	50

Disposable NIBP cuffs



Code (Model)	Usage	Width	Range	Qty
S954A (YP-840P)	Infant	5 cm	8 - 13 cm	20
S954B (YP-841P)	Child	7 cm	13 - 18 cm	20
S954C (YP-842P)	Adult, Small	10 cm	18 - 23 cm	20
S954D (YP-843P)	Adult	13 cm	23 - 33 cm	20
S954E (YP-844P)	Adult, Large	16 cm	33 - 45 cm	20
S954F (YP-845P)	Thigh	19 cm	45 - 55 cm	20



Code (Model)	Usage	Width	Range	Qty
S948A (YP-820P)	Neonate	2 cm	3 - 6 cm	10
S948B (YP-821P)		3 cm	4 - 8 cm	10
S948C (YP-822P)		4 cm	6 - 11 cm	10
S948D (YP-823P)		4.5 cm	7 - 13 cm	10
S948E (YP-824P)		5 cm	8 - 15 cm	10



Code (Model)	Usage	Width	Range	Qty
S945C (YP-810P)	Child	6 cm	8 - 14 cm	20
S945D (YP-811P)		8 cm	13 - 20 cm	20
S946E (YP-812P)	Adult	10 cm	18 - 26 cm	20
S946F (YP-813P)		14 cm	26 - 35 cm	20
S946G (YP-814P)		15 cm	29 - 38 cm	20
S946H (YP-815P)		17 cm	32 - 42 cm	20
S946I (YP-816P)		18 cm	35 - 44 cm	20
S946J (YP-817P)		20 cm	42 - 50 cm	20

YAWARA CUFF2



Yawara(kai)
means
“soft in touch”



Code (Model)	Usage	Width	Range
S951A (YP-710T)	Infant	5 cm	8 - 13 cm
S951B (YP-711T)	Child	7 cm	13 - 18 cm
S951C (YP-712T)	Adult, small	10 cm	18 - 23 cm
S951D (YP-713T)	Adult	13 cm	23 - 33 cm
S951E (YP-714T)	Adult, large	16 cm	33 - 45 cm
S951F (YP-715T)	Thigh	19 cm	45 - 55 cm

Reusable NIBP cuffs



Code (Model)	Usage	Width	Range
S943A (YP-960T)	Infant	5 cm	8 - 13 cm
S943B (YP-961T)	Child	7 cm	12 - 18 cm
S943C (YP-962T)		10 cm	15 - 23 cm
S944B (YP-963T)	Adult	13 cm	21 - 30 cm
S944C (YP-964T)		15 cm	26 - 36 cm
S944D (YP-965T)	Thigh	19 cm	33 - 45 cm

NIBP adapters

Adult/child



Neonate




Code (Model)	Usage	Cord
S902 (YN-901P)	Adult/child	3.5 m
S901 (YN-900P)		1.5 m
S905 (YN-921P)	Neonate	3.5 m
S904 (YN-920P)		1.5 m

Extension cord for S901

Code (Model)	Cord
S903 (YN-990P)	1.5 m

Temperature

Temperature probes



Code (Model)	Type
P242D (409J)	Disk
P241B (402J)	Catheter, child
P240B (401J)	Catheter, adult

Probe covers



Code (Model)	Qty
P249A (-)	10

*Use with P240B

Insulation pads



Code (Model)	Qty
P252 (-)	60

*Use with P242D

Temp connection cord



Code (Model)
K961 (JT-900P)

IBP

IBP connection cords



Code (Model)
K951 (JP-900P)

*Use with MERIT MEDICAL transducer



Code (Model)	Cord
L901 (JP-920P)	3.5 m
L902 (JP-921P)	1.5 m

*Use with Edwards Lifescience transducer



Code (Model)
K957 (JP-960P)

*Use with Biosensors transducer



Code (Model)
K952 (JP-910P)

*5P connector

IBP interface isolation cable

Code (Model)
- (JP-911P)

*Connects an Integra Neuro Sciences Camino® MPM-1 multi-parameter monitor to a BSM-9101K bedside monitor

Multigas (GF-210R/220R)

Water trap



Code (Model)	Qty
V918 (YG-600P)	12

Sampling line



Code (Model)	Qty
V919 (YG-610P)	25

T-Piece




Code (Model)	Type	Dead space	Qty
V920A (YG-620P)	Straight	5 mL	25
V920B (YG-621P)	Elbow	7 mL	25



EEG

Disk electrode, EEG



Code (Model)
H471 (BE-910P)

Collodion electrode, EEG



Code (Model)
H472 (BE-920P)

Electrode lead, EEG



Code (Model)
K461 (BE-930P)

EEG connection cord



Code (Model)
L461 (JE-905P)

EEG disk electrode



Code (Model)
H503B (NE-114A)

EEG collodion electrode



Code (Model)
H527 (NE-136A)

EEG needle electrode



Code (Model)
H537A (NE-224S)

EEG connection cord



Code (Model)
L462 (JE-906P)

EEG Headband / Disposable Electrode for BSM-6000/3000, CSM-1900+AE-918P



Code (Model)	Size
H564B (BE-090AS)	S
H564C (BE-090AM)	M



Code (Model)
H548A (NE-091S7)
*7 pcs × 6 packs/box



Code (Model)
H548B (NE-090S1)
*1 pcs × 6 packs/box

Monitor Flow/Paw, CO, Batteries, Recording Paper, Others

Flow/Paw

Flow adapter



Code (Model)
P913 (TF-120P)

*Use with GF-120P/220P

Flow tube



Code (Model)
P914 (YF-120P)

*Use with GF-120P/220P

Luer cap



Code (Model)
P915 (YS-094P0)

CO

CO connection cord



Code (Model)
K962 (JT-950P)

Batteries



for BSM-3000/6000, WEP-5200

Code (Model)
X075 (SB-671P)



for OLV-2700/OLG-2800

Code (Model)
X072 (YS-085P6)



for PVM-2700

Code (Model)
X076 (SB-201P)



for BSM-1700, OLV-4000, OLG-3800

Code (Model)
X161 (SB-170P)

ECG/BP output cables

K974 (YJ-910P)



Code (Model)
K974 (YJ-910P)

K975 (YJ-920P)



With connector
With connector

Touch pen



for BSM/WEP

Code (Model)
Y075 (-)

Power cord

Code (Model)	
L935	US version, 100 - 110V
L936	European version, 200 - 240V
L944	L-shape, 100 - 127V
L948	L-shape, 200 - 240V

Recording paper

for BSM-2300/4100/5100, WEP-4200, CNS-9701

Code (Model)	Qty	Type	Width × Length
A723 (FQW50-3-100)	10	Z-fold	50 mm × 30 m

for BSM-3000/6000/9100, CSM-1901, CNS-6201/9601, PVM-2700, WEP-5200

Code (Model)	Qty	Type	Width × Length
A721 (FQW50-2-100)	10	Z-fold	50 mm × 20 m

for OPV-1500






Code (Model)	Qty	Type	Width × Length
A810 (QW60-2)	10	Roll	60 mm × 20 m

Thermal head cleaning pens










Code (Model)	Qty
Y011 (-)	5

Monitor Accessory Set

			For BSM-6000/3500/3700/1700, CSM-1900 + NK SpO ₂ probe				For BSM-6000/3500/3700/1700, CSM-1900 + Masimo / Nellcor SpO ₂ probe	
Order code			Y212A	Y212B	Y213A	Y213B	Y212C	Y212D
Model		Qty	YO-60AY1	YO-60AY2	YO-60AY1	YO-60AY2	YO-60AY3	YO-60AY4
Electrode lead 3 electrode, clip, 0.8 m		1	IEC type (Red, Yellow, Green) K911 (BR-903P)		AHA type (White, Black, Red) K911A (BR-903PA)		IEC type (Red, Yellow, Green) K911 (BR-903P)	
ECG connection cord, 3 m		1	IEC type K922 (JC-906P)		AHA type K922A (JC-906PA)		IEC type K922 (JC-906P)	
SpO ₂ connection cord, 2.5 m		1	K931 (JL-900P)					
Air hose for NIBP, 3.5 m		1	S902 (YN-901P)					
NIBP cuff		1	13 cm width 23-33 cm range S951D (YP-713T)	16 cm width 33-45 cm range S951E (YP-714T)	13 cm width 23-33 cm range S951D (YP-713T)	16 cm width 33-45 cm range S951E (YP-714T)	13 cm width 23-33 cm range S951D (YP-713T)	16 cm width 33-45 cm range S951E (YP-714T)

Power cord is not included in the set. Order code: L935 (H type, 100-110V), L936 (N type, 200-240V)
The model of cuff might be different depending on the shipment time.

			For PVM-2700			
Order code			YS-101P0	YS-101P1	YS-101P2	YS-101P3
		Qty				
Electrode lead 3 electrode, clip, 0.8 m		1	IEC type (Red, Yellow, Green) K911 (BR-903P)		AHA type (White, Black, Red) K911A (BR-903PA)	
ECG connection cord, 3 m		1	IEC type K922 (JC-906P)		AHA type K922A (JC-906PA)	
Disposable electrode general use, 35mm dia		150	150 pcs/box G207 (L-150X)			
SpO ₂ connection cord, 2.5 m		1	K931 (JL-900P)			
Air hose for NIBP, 3.5 m		1	S902 (YN-901P)			
NIBP cuff		1	13 cm width 21-30 cm range S944B (YP-963T)	15 cm width 26-36 cm range S944C (YP-964T)	13 cm width 21-30 cm range S944B (YP-963T)	15 cm width 26-36 cm range S944C (YP-964T)
Battery		1	X076 (SB-201P)			

Power cord is not included in the set. Order code: L935 (H type, 100-110V), L936 (N type, 200-240V)